

# 



pH ORP Ion Conductivity

Resistivity Total Dissolved Solids Salinity

Benchtop Water Quality Instruments
Colour Touchscreen Meters





www.horiba-laqua.com



# Benchtop Water Quality Instruments

Colour Touchscreen Meters



#### 2003

F-50 (desktop) The world's first pH meter with colour LCD display. Navigation panel guides operators on how to use the meter as well as resolve errors.



D-50 (portable) Waterproof, IP67rated housing and multi-parameter.

#### 2011



LAQUA Benchtop Water Quality Instruments

LAQUAtwin Pocket Water Quality Meters



LAQUA Handheld Water Quality Instruments

HORIBA

#### 1993

F-20 (benchtop) The world's first wireless pH meter. Large graphical display gives user instructions on screen.

#### 1990

B-111 (Pen type) The pen type sensor allows small samples to be tested.



#### 1987

C-1 (card) Development of the world's first flat sensor.



#### 1980

Model F-80 (benchtop) The world's first instrument capable of measuring pH at 0.001 resolution includes an integral computer with automatic calibration and a self-diagnostic function.

L-7 (integrated) Introduction of a small, handheld pH meter with integrated electrode.

# History of the HORIBA pH Meter

#### The humble beginning of HORIBA...

In 1950, Dr. Masao Horiba pioneered and launched Asia's first pH meter in Kyoto, Japan. Since then, HORIBA has been introducing several of the world's firsts such as the first 0.001 resolution pH meter, the first flat sensor featured in the Cardy, the first wireless pH meter, the first colour LCD display, etc.



Model F-7AD (benchtop) Incorporating an industry-first LCD display, the combination of a glass electrode, a reference electrode and a temperature-compensating electrode, makes testing easier.



#### 1964

M-5 (benchtop) conversion from vacuum tube to semiconductor allows miniaturization and development of fast response meter



### 1950

**HORIBA** introduces Japan's first glass electrode pH meter.



- Large touch screen color graphic LCD—5.7 inches (115.2 x 86.4 mm)
- Chemical-resistant, 2mm thick super white glass panel with protection cover
- Easy to clean and elegant round body
- GLP / GMP compliant
- Switchable display—digital, graph, and analog

- Effortless single-touch operations—tap, flick, and drag
- 2-Channel display and simultaneous measurements for F-73 and F-74 models
- Small footprint—170 (W) x 174 (D) x 73 (H) mm
- Data acquisition software in mini USB is included
- 21 CFR Part 11 software complies with U.S. FDA's system requirements for electronic records and signatures (optional)



Protection Cover



Data Acquisition Software



21 CFR Part 11 Software



# Intuitive Touch-Control Operation



6 types of international standard plugs included (US, UK, EU, Australia / New Zealand, Korea and China)

#### Meter Connections

# Data Management

#### Data Key




 Data key shows settings that allow users to search, view, delete, and copy data from meter to USB flash drive

#### Sample ID



 Meter internal memory stores up to 2000 data with sample ID for easy reference

#### Data Storage



- Data can be stored simultaneously on both meter and USB flash drive (if inserted)
- Calibration and measurement data are logged automatically at set time interval

#### Data Search




• Data search by date, parameter, or sample ID



- Data output via USB to PC / USB flash drive or via RS232C to PC / printer (cables sold separately)
- Analog output adjustment—voltage output can be acquired from digital multimeter or recorder connected to the analog output connector

#### **Custom Printout**



- Auto or manual printing of calibration and measurement values for record keeping
- Printout contents can be customized based on user preference or GMP/GLP requirements—date and time, operator, electrode and meter information, electrode status, and calibration data

#### Meter Security



- Password setting for security
- Up to 25 administrators or operators can be registered



### Calibration Support Function

Enjoy hassle-free calibration with on screen support. The meter will walk you through the steps of calibration.

- Auto Buffer Recognition
- Auto Calibration Function







# Reading Stability Check

- · Perform proper calibration with stable readings
- Determine the stability of reading at a glance in either digital or graph display during pH and ion calibration
- Stability value is a deviation between the maximum and minimum readings in the last 10 seconds





#### Electrode Status

- Electrode condition and results such as calibrated values, offset, acid and alkaline slopes, are shown at the end of calibration
- Programmable calibration reminder and alarm for measured values exceeding set limits
- Temperature indicator appears when a temperature probe or electrode with integrated temperature sensor is connected to the meter
- Temperature sensor calibration function
- Electrode model, either selected from preset list or entered manually, and lot or MFG no. (entered manually) are included in stored data and printouts

# Inspection Function

Easy navigation for meter and electrode inspections using a simulator. Various industrial standards (JIS, USP, EP, JP, CP) are also supported.

Convenient for IQ / OQ / PQ validation





# NAVIGATION

### **Troubleshooting Function**

On-screen support for resolving a problem that occurs during calibration or sample measurements. A user's guide is incorporated in the software to assist with any operational difficulties.




4

HCP

EP.


CP



# **Application Functions**

Various industry standard methods are supported by the instrument. Conductivity measurement for several pharmaceutical pure water guidelines and ion standard addition methods are incorporated in the meter.






>

>

2

COND pharmacopoela mode



standard method. Just wait for the result.

(C) MEAS

2016/05/18 17:44



#### На

- 5 pH buffer groups
  - o USA (1.68, 4.01, 7.00, 10.01, 12.45)
  - o NIST (1.68, 4.01, 6.86, 9.18, 12.45)
  - o NIST2 (1.68, 4.01, 6.86, 10.01, 12.45)
  - o China (1.68, 4.01, 6.86, 9.18, 12.46)
  - Custom (any pH buffers)

- SET CHI 0.01 pH 0.001 pH TEMP conversion Auto Atlanta, upper Italia Alarm, Lower Limit OFF. W. Electrode:model Customize | -
- SET CHI D MIST USA NIST2 China CUSTON

- Up to 5 calibration points
- 0.01 and 0.001 pH Resolutions
- Auto setting allows the meter to toggle between 0.01 and 0.001 resolution depending on the stability of the reading
- Auto calibration / Auto buffer recognition

Display absolute potential and relative potential





# ADVANCED

CH1

ORP

8 0 m

#### **ORP**

Capable of 1-point calibration

# lon

- Make your own calibration curve with maximum of 5 points or perform standard addition techniques
- Programmed with standard addition methods—known addition and sample addition (single and double are available for both methods)
- Measurement units µg/L, mg/L, g/L, mmol/L, mol/L

# ORP 8 0 M




3

>

>

>

#### 2016/06/27 11:21 ION standard addition mode Known addition method > Sample addition method > TEMP setting ATC |



# Conductivity

- Automatic / manual calibration up to 4 points
- Adjustable temperature coefficient and reference temperature for temperature compensated readings
- Selectable cell constants 0.1, 1.0, 10.0
- Auto ranging S/cm and S/m units, fix mS/cm unit
- Support conductivity standard methods for pharmaceutical water—USP, EP, JP and CP

# Total Dissolved Solids (TDS)

- Programmed with 4 predetermined TDS curves for accurate measurement—Linear, EN27888, 442, and NaCl
- Select the TDS curve suitable for your application
- Calibration only in conductivity mode is required

#### **TDS Calibration Curves**

MEAS

CH2 ►

COND >


	Application	Key chemical species	TDS selection
	Aquaculture, pickling	NaCl	NaCl
	Boiler water, HVAC	Na <sub>2</sub> SO <sub>4</sub> , NaHCO <sub>3</sub> , NaCl	442 (Myron)
	Environmental	EN standard for environmental water	EN 27888
	General application	Not known	KCI (linear factor) Default: 0.5 Selectable: 0.4 to 1.0

### Salinity

- Programmed with 2 predetermined salinity curves—NaCl and seawater
- Salinity value is calculated based on measured conductivity value
- 1-point calibration using standard solution
- Measurement units—percentage (%) and parts per thousand (ppt)







#### Auto Stable / Auto Hold

- In measurement mode, the meter displays live readings continuously
- Activate auto hold by tapping START
- Auto hold settings—Exact, Normal, Brief, Time, Customize, and Manual





# FEATURES

### Auto Log Data

 Log data automatically by setting time interval from 1 to 999 seconds





# Multi-Language

 Choose a language that you are familiar with—English, Japanese, Chinese, Korean, and Vietnamese




# Screen Settings

- Set stylish theme on your meter screen—Standard, Cool, Monotone, and Kyoto
- Power saving mode—turns off the backlight to save power

# Sound Setting


 Play a click sound every time you tap a key







Accessories				
Code	Part No.	Description		
LAQUA-SW-21CFR11	3200707161	21 CFR Part 11 Software includes CD with PIN code, USB cable, and manual		
Printer Printer	3014030147 (230v) 3014030146 (120v)	Printer (for GLP/GMP compliance) Cable sold separately, Plain paper		
cable	3014030148	Printer cable (1.5 m)		
Ink Printer	3014030149	Printer paper (20 rolls)		
ribbon paper	3014030150	Ink ribbon (5 pcs/set)		
Universal AC adapter	3200647413	Multi-Voltage (100-240V) with 6 plugs, (US, UK, EU, ANZ, Korea and China) 1.8 m cable		
686,	3014028368	Digital simulator X-51 (pH, mV, Ion, DO, temperature simulator)		
X-51 X-52	3014028370	Digital simulator X-52 (Conductivity, temperature simulator)		
	3200382462	LCD protection sheet (2 pcs/pack)		
LCD Protection protection cover sheet	3200382441	Protection cover (Protects the meter for F-70, DS-70, 1000 series)		
	3200373941	USB cable (to connect meter and PC.)		
9	3014030152	Analog cable (Analog (alarm) output cable)		
USB Serial cable	3014030151	Serial cable (to connect meter and PC (Serial, 9 pins))		
FA-70S	3200382557	Adjustable, free-standing electrode stand (Height: 384 mm) image on the right		
FA-70L	3200382560	Long, free-standing electrode stand (Height: 450-650mm) image on the right		
- 1215	3200373991	Arm for electrode stand FA-70A, FA-70S, & FA-70L		
	3200373961	Electrode holders, 2pcs (for mounting electrode with round cap on electrode stand arm)		
	3200382477	Electrode protection caps, 3pcs (for 9615S-10D, 9618S-10D, 9681S-10D pH electrode)		
	3200043508	Electrode protection caps, 5pcs (for 9621-10D, 9625-10D, 9630-10D, 9631-10D, 9632-10D, 6367-10D, 6377-10D, 6252-10D, 6261-10C, 1066A-10C, 1076-10C, 2060-10T, 9300-10D, 9382-10D, 3552-10D pH electrode)		
	3200382482	Electrode protection cap for long electrode (for 9680S-10D, 9480-10C pH Electrode)		



With over 60 years of engineering excellence, HORIBA's diverse range of water quality analyzers and electrodes are ideal for everyday laboratory needs through to the most demanding of applications. Visit our website for a wealth of useful information and water quality measurement tips to help you obtain the best results in your work.







#### **Benchtop Meters**

Developed using extensive feedback from users, our new LAQUA meters deliver the best solution for water quality analysis. Our LAQUA website features an online 'Selection Guide' to enable you to find the perfect LAQUA meter and electrode for your need.



#### Handheld Meters

In the lab, in the field or anywhere you need it. LAQUA Handheld meters are designed for use with one hand and with an IP67 waterproof rating and shock-resistant casing. Meters can be used for long periods, even in dark places, making it ideal for field measurements in rivers and lakes.



#### **Electrodes**

Various electrodes to match any application. A wide range of products for both benchtop and portable systems are available, including easy and reliable standard models, applicationfocused models for small samples or large containers, and special electrodes for specific sample characteristics.



#### **Pocket Meters**

Analyzing water quality is simplified when using our LAQUAtwin range of meters. Designed to produce accurate and reliable results. Anyone, anywhere, at any time can measure samples easily with a LAQUAtwin meter. See just how good they are at our website.





#### **Application Notes**

LAQUAtwin pocket meters offer quick and convenient alternative to analyze important parameters with high accuracy. Several application notes are available at (http://goo.gl/znwE6j) detailing the use of LAQUAtwin and the results achieved for the respective applications. Additional application notes will be added when available.













- The contents of this catalog are subject to change without prior notice, and without any subsequent liability to this company.
- The color of the actual products may differ from the color pictured in this catalog due to printing limitations
- It is strictly forbidden to copy the content of this catalog in part or in full.
- All brand names, product names and service names in this catalog are trademarks or registered trademarks of their respective companies.
- Windows is a registered trademark of Microsoft Corporation in the United States and other countries

#### HORIBA Instruments (Singapore) Pte. Ltd.

83 Science Park Drive, #02-02A, The Curie, Singapore 118258 Phone: 65 6908-9660 Fax: 65 6745-8155

e-mail: laqua@horiba.com www.horiba-laqua.com

#### **HORIBA UK Limited**

Kyoto Close, Moulton Park, Northampton NN3 6FL Phone: 44 (0) 1604 542567 Fax: 44 (0) 1604 542699

e-mail: waterquality@horiba.com www.horiba.com/uk



Brochure HBTC-05-2017B

