



Free Download

<https://updates.waterlinkconnect.com/stable/WaterLinkConnect.exe>

WaterLink[®] Connect 2 Software for Transferring your Water Testing Results

Connect your meter. Collect your results.

WaterLink Connect v2.1.0.0

Meter Controls

Location
Type here to add or search...

Sample
Salt

Reagent
203 Series Cl/Br With Phosphate Disk

Start Test

Retrieve Log

Results

Filter by Time

From: 08-05-2000 02:18:18 AM
To: 03-15-2018 12:37:07 PM

Filter by Test Name

Alkalinity
 Combined Cl
 Copper
 Free Cl
 Hardness

Filter by Location

None

Filter by Meter

Spin Touch

Timestamp	Name	Value	Units	Meter	Location
3/15/2018 12:37 PM	Free Cl	2.26	ppm	Spin Touch	None
3/15/2018 12:37 PM	Phosphate	2521	ppb	Spin Touch	None
3/15/2018 12:37 PM	Iron	1.3	ppm	Spin Touch	None
3/15/2018 12:37 PM	Copper	4	ppm	Spin Touch	None
3/15/2018 12:37 PM	Stabilizer	39	ppm	Spin Touch	None
3/15/2018 12:37 PM	Hardness	0	ppm	Spin Touch	None
3/15/2018 12:37 PM	Alkalinity	97	ppm	Spin Touch	None
3/15/2018 12:37 PM	pH	7.4		Spin Touch	None
3/15/2018 12:37 PM	Combined Cl	0	ppm	Spin Touch	None
3/15/2018 12:37 PM	Total Cl	2.26	ppm	Spin Touch	None
3/14/2018 12:51 PM	Copper	0	ppm	Spin Touch	None
3/14/2018 12:51 PM	Stabilizer	5	ppm	Spin Touch	None
3/14/2018 12:51 PM	Hardness	0	ppm	Spin Touch	None
3/14/2018 12:51 PM	Alkalinity	0	ppm	Spin Touch	None
3/14/2018 12:51 PM	Phosphate	0	ppb	Spin Touch	None
3/14/2018 12:51 PM	Combined Cl	0	ppm	Spin Touch	None
3/14/2018 12:51 PM	Total Cl	0	ppm	Spin Touch	None
3/14/2018 12:51 PM	Free Cl	0	ppm	Spin Touch	None
3/14/2018 12:51 PM	Iron	0	ppm	Spin Touch	None
3/14/2018 12:51 PM	pH	6.3		Spin Touch	None
3/14/2018 8:50 AM	Iron	1.3	ppm	Spin Touch	None
3/14/2018 8:50 AM	Phosphate	2520	ppb	Spin Touch	None

WaterLink® Connect 2 Software

Supporting some of LaMotte Company's most advanced meters, **WaterLink® Connect 2** is a free Windows® application for collecting test results over a USB connection and storing them in a local database. Assign test results to a testing location, or filter your results and export them to a CSV format. Additionally, WaterLink® Connect acts as a bridge between your meter and detailed water analysis solutions like DataMate® Web (datamateweb.php) and WaterLink® Solutions™ (waterlinksolutions.php) .

- Collect numeric test results from supported meters in real time testing, or as a bulk transfer*
- Filter results by date, test factor name, meter type, and assigned location
- Export results to CSV

WaterLink® Connect 2 supports USB connection to the following LaMotte meters:

- WaterLink® Spin
- WaterLink® Spin Touch® (<http://www.waterlinkspintouch.com/>)
- DC1500 (<http://www.lamotte.com/en/industrial/instrumentation/3240.html>)
- SMART3 (<http://www.lamotte.com/en/water-wastewater/instrumentation/1910.html/>)

*Result collection methods vary based on meter capabilities.

Free Tech Support Information

ADDRESS

802 Washington Ave.
Chestertown, Maryland 21620 USA

PHONE

USA: (800) 344-3100 (Option 3)

EMAIL

softwaresupport@lamotte.com
(mailto:softwaresupport@lamotte.com)

OFFICE HOURS

Mon - Fri / 9:00 AM - 5:00 PM (EST)

LaMotte Software Center

WaterLink DataMate Web Software (<http://softwarecenter.lamotte.com/datamateweb.php>)

WaterLink Solutions Software (<http://softwarecenter.lamotte.com/waterlinksolutions.php>)

WaterLink Connect 1 Software (<http://softwarecenter/waterlinkconnect1.php>)

WaterLink Connect 2 Software (<http://softwarecenter.lamotte.com/waterlinkconnect2.php>)

Support (<http://www.lamotte.com/en/remotesupport>)

© Copyright 2018 LaMotte Company. All Rights Reserved. / Privacy Policy

(<http://softwarecenter.lamotte.com/privacy.php>) / Terms and Conditions (<http://softwarecenter.lamotte.com/tos.php>)